



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 75 x 13 mm, with paper label**

Page 1



Product description

Order number	04.1919.100
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 3 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-28

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 75 x 13 mm, with paper label**

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 3 ml
Diameter 13 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap blue
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitro diagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of box 171 mm
Width of box 116 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 75 x 13 mm, with paper label**

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0045 kg
Weight of case	3.01 kg
EAN of inner box	4038917266970
EAN case	4038917266963

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

