



**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with plastic label**

**Page 1**



**Product description**

Order number	04.1925
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 4 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Serum Gel CAT, 4 ml, cap brown,  
(LxØ): 75 x 13 mm, with plastic label**

**Page 2**

Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	4 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with plastic label**

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0054 kg
Weight of case	3.47 kg
EAN of inner box	4038917248310
EAN case	4038917248303

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

