



Product Specification
S-Monovette® neutral Z, 4.9 ml, cap white, (LxØ):
90 x 13 mm, with paper label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 04.1926.001 |
| Product description | S-Monovette® neutral Z, without preparation, 4.9 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: black, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|---------------------|
| Type of collection | venous |
| Type of preparation | without preparation |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |
| Liquid-tight | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® neutral Z, 4.9 ml, cap white, (LxØ):
90 x 13 mm, with paper label

| | |
|-------------------|---------------|
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 150 mm |
| Case volume | 0.0208 cbm |
| Weight of product | 0.0048 kg |
| Weight of case | 3.18 kg |
| EAN of inner box | 4038917134422 |
| EAN case | 4038917025683 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

