



### Product description

Order number	04.1929
Product description	S-Monovette® Lithium heparin LH, LH, preparation: Lithium heparin, 2.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 2.7 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label**

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0045 kg
Weight of case	3 kg
EAN of inner box	4038917360678
EAN case	4038917360661

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

