



## Product Specification

**S-Monovette® Fluoride/EDTA FE, 4.9 ml, cap grey,  
(LxØ): 90 x 13 mm, with paper label**

**Page 1**

### Product description

Order number	04.1933.100
Product description	S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 4.9 ml, membrane screw cap, cap grey, colour code ISO, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: grey, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Fluoride + EDTA
Preparation concentration	1.0 mg/ml blood + 1.2 mg/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	grey
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	4.9 ml
Diameter	13 mm
Length including cap	106 mm
Length excluding cap	90 mm
Length including cap and plunger	127 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	grey
Membrane material	rubber
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Fluoride/EDTA FE, 4.9 ml, cap grey,  
(LxØ): 90 x 13 mm, with paper label**

**Page 2**

Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	36000
Depth of box	171 mm
Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0049 kg
Weight of case	3.22 kg
EAN of inner box	4038917366441
EAN case	4038917366434

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**