



Product description

Order number	04.1935
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 4.9 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 90 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum Gel CAT, 4.9 ml, cap brown,
(LxØ): 90 x 13 mm, with plastic label

Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0057 kg
Weight of case	3.62 kg
EAN of inner box	4038917088077
EAN case	4038917006217

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

