



**Product Specification**  
**S-Monovette® Lithium heparin LH, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label**

**Page 1**



**Product description**

Order number	04.1936.001
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 4.9 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: orange, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label**

**Page 2**

Liquid-tight Connection                      yes  
membrane screw cap

**Size**

Sample volume                                  4.9 ml  
Diameter    13 mm  
Length including cap                          106 mm  
Length excluding cap                          90 mm  
Length including cap and plunger           127 mm

**Material & colours**

Colour of product                              transparent  
Closure material                                High Density Polyethylene (HD-PE)  
Colour of cap                                    orange  
Membrane material                            rubber  
Tube material                                    Polypropylene (PP)  
Piston material                                  High Density Polyethylene (HD-PE)  
Piston plunger material                      Polystyrene (PS)

**Purity & certification**

Satisfies the requirement                    IATA, ADR  
Product category                                In vitrodiagnostic device, CE  
CE certified                                      CE  
Purity standard                                  sterile  
Sterilisation                                      Electron irradiation  
Batched    yes

**Packaging**

Minimum order qty.                            500  
Type of smallest subpackaging              case  
Piece(s) / inner box                            50  
Piece (s) / outer case                          500  
Piece(s) / pallet                                36000  
Depth of box                                    171 mm  
Width of box                                      130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label**

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0054 kg
Weight of case	3.45 kg
EAN of inner box	4038917134743
EAN case	4038917025911

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

