



**Product Specification**  
**S-Monovette® Serum CAT, 2.7 ml, cap red, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 1**



**Product description**

Order number	04.1943.100
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 2.7 ml, membrane screw cap, cap red, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: red, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







**Product Specification**  
**S-Monovette® Serum CAT, 2.7 ml, cap red, (LxØ):**  
**75 x 13 mm, with paper label**

Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0045 kg
Weight of case	3.02 kg
EAN of inner box	4038917266857
EAN case	4038917266840

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

