




Product Specification

**S-Monovette® EDTA K3E, 1.8 ml, cap red, (LxØ):
65 x 13 mm, with plastic label pre-barcoded,
Pre-barcode with 8-digit unique number range
and 3-digit prefix**

Page 1



Product description

Order number	04.1951.088
Product description	S-Monovette® EDTA K3E, K3E, preparation: K3 EDTA, 1.8 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 013, label/print: white/red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® EDTA K3E, 1.8 ml, cap red, (LxØ):
65 x 13 mm, with plastic label pre-barcoded,
Pre-barcode with 8-digit unique number range
and 3-digit prefix**

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 1.8 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 81 mm
Length excluding cap 65 mm
Length including cap and plunger 102 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 20000
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 42000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® EDTA K3E, 1.8 ml, cap red, (LxØ):
65 x 13 mm, with plastic label pre-barcoded,
Pre-barcode with 8-digit unique number range
and 3-digit prefix**

Page 3

Depth of box	171 mm
Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0039 kg
Weight of case	2.63 kg
EAN of inner box	4038917949385
EAN case	4038917949392

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT