



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml,  
cap orange, (LxØ): 75 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1952.001
Product description	S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 2.7 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel+
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	7 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-28

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml,  
cap orange, (LxØ): 75 x 13 mm, with paper label**

**Page 2**

Liquid-tight yes  
Connection membrane screw cap

### Size

Sample volume 2.7 ml  
Diameter 13 mm  
Length including cap 91 mm  
Length excluding cap 75 mm  
Length including cap and plunger 112 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap orange  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitro diagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 39000  
Depth of box 171 mm  
Width of box 116 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



Product Specification

**S-Monovette® Lithium heparin gel+ LH, 2.7 ml,  
cap orange, (LxØ): 75 x 13 mm, with paper label**

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0051 kg
Weight of case	3.31 kg
EAN of inner box	4038917523035
EAN case	4038917523028

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

