



Product Specification
**S-Monovette® Serum CAT, 1.6 ml, cap red, (LxØ):
75 x 13 mm, with plastic label**

Page 1



Product description

Order number	04.2911.200
Product description	S-Monovette® Serum CAT, CAT, preparation: Clotting activator, 1.6 ml, membrane screw cap, cap red, colour code ISO, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	transparent/red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Serum CAT, 1.6 ml, cap red, (LxØ):
75 x 13 mm, with plastic label

Page 2

Connection membrane screw cap

Size

Sample volume 1.6 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of box 171 mm
Width of box 116 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 1.6 ml, cap red, (LxØ):
75 x 13 mm, with plastic label

Page 3

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0049 kg
Weight of case	3.18 kg
EAN of inner box	4038917983327
EAN case	4038917983334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT