



Product description

Order number	04.2913
Product description	S-Monovette® Lithium heparin LH, LH, preparation: Lithium heparin, 1.6 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 1.6 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label

Page 2

Liquid-tight Connection yes
membrane screw cap

Size

Sample volume 1.6 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap orange
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging case
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 1.6 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0051 kg
Weight of case	3.28 kg
EAN of inner box	4038917983389
EAN case	4038917983396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

