



Product description

Order number	04.2920.781
Product description	, membrane screw cap, (LxØ) without cap: 65 x 13 mm, with plastic label, sterile

Product characteristics

Type of collection	venous
Preparation concentration	1.6 mg/ml blood
Label/ Print	with plastic label
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

Size

Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	65 mm
Length including cap and plunger	102 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Purity standard	sterile
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Piece(s) / inner box	50
Piece (s) / outer case	500
Depth of box	170 mm
Width of box	115 mm
Height of box	72 mm
Depth of case	380 mm
Width of case	350 mm
Height of case	120 mm
Case volume	0.016 cbm
EAN of inner box	4038917984010
EAN case	4038917984027

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

