



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,  
cap blue, (LxØ): 66 x 11 mm, with paper label**

**Page 1**



### Product description

Order number	05.1165.100
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 3 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,  
cap blue, (LxØ): 66 x 11 mm, with paper label**

**Page 2**

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

### Size

Sample volume	3 ml
Outer diameter of cap	13 mm
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	56000
Depth of box	145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,  
cap blue, (LxØ): 66 x 11 mm, with paper label**

**Page 3**

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0041 kg
Weight of case	2.61 kg
EAN of inner box	4038917126113
EAN case	4038917020442

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**