



Product Specification
S-Monovette® EDTA K2E, 1.2 ml, cap violet,
(LxØ): 66 x 8 mm, with plastic label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 06.1660.100 |
| Product description | S-Monovette® EDTA K2E, preparation: K2 EDTA, 1.2 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/violet, 50 piece(s)/case, sterile |

Product characteristics

| | |
|---------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | K2 EDTA |
| Preparation concentration | 1.6 mg/ml blood |
| Colour code | ISO |
| Label/ Print | with plastic label |
| Colour of print/label | white/violet |
| Graduation | yes |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Liquid-tight | yes |
| Connection | membrane screw cap |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® EDTA K2E, 1.2 ml, cap violet,
(LxØ): 66 x 8 mm, with plastic label

Page 2

Size

| | |
|----------------------------------|--------|
| Sample volume | 1.2 ml |
| Diameter | 8 mm |
| Length including cap | 82 mm |
| Length excluding cap | 66 mm |
| Length including cap and plunger | 103 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | violet |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | rubber |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 56000 |
| Depth of box | 145 mm |
| Width of box | 108 mm |
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® EDTA K2E, 1.2 ml, cap violet,
(LxØ): 66 x 8 mm, with plastic label

| | |
|-------------------|---------------|
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.0033 kg |
| Weight of case | 2.18 kg |
| EAN of inner box | 4038917322737 |
| EAN case | 4038917322720 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

