



Product description

Order number	02.1063.014
Product description	S-Monovette® Serum, CAT, preparation: Clotting activator, 9 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 92 x 16 mm, with paper label, label / print: white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with paper label
Colour of print/label	white
Graduation	fill mark
Cap	membrane screw cap
Type of cap	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20 °C
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® Serum, 9 ml, cap white, (LxØ): 92 x 16 mm, with paper label

Page 2

Connection membrane screw cap

Size

Sample volume 9 ml
Diameter 16 mm
Length including cap 108 mm
Length excluding cap 92 mm
Length including cap and plunger 129 mm

Material & colours

Colour of product transparent
Cap material High Density Polyethylene (HD-PE)
Colour of cap white
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of inner packaging case
Piece(s) / inner box 50
Piece(s) / case 500
Piece(s) / pallet 24000
Depth of case 648 mm
Width of case 303 mm
Height of case 150 mm
Case volume 0,0295 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Serum, 9 ml, cap white, (LxØ): 92 x 16 mm, with paper label

Page 3

Weight of product	0,0088 kg
Weight of case	5,35 kg
EAN of inner box	4038917124966
EAN case	4038917019576

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT