



Product Specification
Luer Monovette® neutral, 9 ml, cap natural, (LxØ):
92 x 16 mm, with paper label

Page 1



Product description

Order number	02.226.160
Product description	Luer Monovette® neutral, without preparation, 9 ml, cap natural, (LxØ) without cap: 92 x 16 mm, with paper label, label / print: black, animal venous blood collection, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	black
Cap	Luer screw cap
Type of cap	screw cap
Application	animal venous blood collection
Recommended centrifugation	2000 x g
Centrifugation time	10 Min
Connection	Luer

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Luer Monovette® neutral, 9 ml, cap natural, (LxØ):
92 x 16 mm, with paper label

Size

Sample volume	9 ml
Diameter	16 mm
Height of product	92 mm
Length of product	92 mm
Length including cap	117 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	100
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	20000
Case volume	0,0382 cbm
Weight of product	0,0101 kg
Weight of case	5,3 kg
EAN of inner box	4038917104326
EAN case	4038917032254

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

