




**Product Specification**  
**S-Monovette® Serum CAT, 5.5 ml, cap white,**  
**(LxØ): 75 x 15 mm, with paper label**

**Page 1**



**Product description**

Order number	03.1397.001
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 5.5 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







**Product Specification**  
**S-Monovette® Serum CAT, 5.5 ml, cap white,**  
**(LxØ): 75 x 15 mm, with paper label**

**Page 3**

Height of box	82 mm
Depth of case	400 mm
Width of case	129 mm
Height of case	440 mm
Case volume	0.0227 cbm
Weight of product	0.007 kg
Weight of case	4.28 kg
EAN of inner box	4038917128438
EAN case	4038917021722

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**