



Product Specification
S-Monovette® Serum, 5.5 ml, cap white, (LxØ): 75 x 15 mm, with plastic label

Page 1



Product description

Order number	03.1397
Product description	S-Monovette® Serum, preparation: Clotting activator, 5.5 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 75 x 15 mm, with plastic label, label/print: transparent/white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**S-Monovette® Serum, 5.5 ml, cap white, (LxØ): 75
x 15 mm, with plastic label**

Page 3

Depth of case	400 mm
Width of case	129 mm
Height of case	440 mm
Case volume	0.0227 cbm
Weight of product	0.0069 kg
Weight of case	4.26 kg
EAN of inner box	4038917086479
EAN case	4038917005036

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT