



**Product Specification**  
**S-Monovette® K2 EDTA, 2.7 ml, cap violet, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1915.100
Product description	S-Monovette® K2 EDTA, preparation: K2 EDTA, 2.7 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: violet, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	K2 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	violet
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-12-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® K2 EDTA, 2.7 ml, cap violet, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 2**

### Size

Sample volume	2.7 ml
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	39000
Depth of box	171 mm
Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-12-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® K2 EDTA, 2.7 ml, cap violet, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 3**

Case volume	0.0188 cbm
Weight of product	0.0042 kg
Weight of case	2.86 kg
EAN of inner box	4038917382250
EAN case	4038917382243

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-12-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**