



**Product Specification**  
**S-Monovette® Serum Gel, 4 ml, cap brown, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1925.001
Product description	S-Monovette® Serum Gel, preparation: clotting activator/gel, 4 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	clotting activator/gel
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-10-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Serum Gel, 4 ml, cap brown, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 2**

Connection membrane screw cap

**Size**

Sample volume 4 ml  
Diameter 13 mm  
Length including cap 91 mm  
Length excluding cap 75 mm  
Length including cap and plunger 112 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap brown  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In-vitro diagnostic, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 20000  
Type of inner packaging case  
Piece(s) / inner box 50  
Piece(s) / case 500  
Piece(s) / pallet 39000  
Depth of box 171 mm  
Width of box 116 mm  
Height of box 73 mm  
Depth of case 388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-10-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Serum Gel, 4 ml, cap brown, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 3**

Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0055 kg
Weight of case	3.49 kg
EAN of inner box	4038917359443
EAN case	4038917359436

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-10-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**