



Product Specification

S-Monovette® Fluoride/EDTA FE, 4.9 ml, cap yellow, (LxØ): 90 x 13 mm, with paper label

Page 1

Product description

| | |
|---------------------|--|
| Order number | 04.1933.001 |
| Product description | S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 4.9 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: yellow, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|-----------------------------------|
| Type of collection | venous |
| Type of preparation | Fluoride + EDTA |
| Preparation concentration | 1.0 mg/ml blood + 1.2 mg/ml blood |
| Colour code | EU |
| Label/ Print | with paper label |
| Colour of print/label | yellow |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2000 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |
| Liquid-tight | yes |
| Connection | membrane screw cap |

Size

| | |
|----------------------------------|--------|
| Sample volume | 4.9 ml |
| Diameter | 13 mm |
| Length including cap | 106 mm |
| Length excluding cap | 90 mm |
| Length including cap and plunger | 127 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Fluoride/EDTA FE, 4.9 ml, cap yellow, (LxØ): 90 x 13 mm, with paper label

Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 20000
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 36000
Depth of box 171 mm
Width of box 130 mm
Height of box 73 mm
Depth of case 388 mm
Width of case 358 mm
Height of case 150 mm
Case volume 0.0208 cbm
Weight of product 0.0049 kg
Weight of case 3.22 kg
EAN of inner box 4038917245968
EAN case 4038917245951

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

