



Product Specification
S-Monovette® DNA Exact, 2.7 ml, cap blue, (LxØ):
75 x 13 mm, with paper label

Page 1



Product description

Order number	04.1948.001
Product description	S-Monovette® DNA Exact, preparation: EDTA based gDNA stabilizer, 2.7 ml, membrane screw cap, cap blue, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: blue, venous blood collection, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	EDTA based gDNA stabilizer
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Design	mit validierten Zeit- und Temperaturparametern
Application	venous blood collection
Example of use	geeignet für die Stabilisierung von genomischer DNA

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® DNA Exact, 2.7 ml, cap blue, (LxØ):
75 x 13 mm, with paper label

Page 2

Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	2.7 ml
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	50
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	39000
Depth of box	171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® DNA Exact, 2.7 ml, cap blue, (LxØ):
75 x 13 mm, with paper label

Page 3

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0042 kg
Weight of case	2.86 kg
EAN of inner box	4038917370745
EAN case	4038917370738

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT