



Product Specification

S-Monovette® Lithium heparin gel+, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label

Page 1



Product description

Order number	04.1954.001
Product description	S-Monovette® Lithium heparin gel+, preparation: Gel+/heparin, 4.9 ml, membrane screw cap, cap orange, colour code ISO, (LxØ) without cap: 90 x 13 mm, with paper label, label / print: brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Gel+/heparin
Preparation concentration	25 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Type of cap	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20 °C
Centrifugation time	7 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette[®] Lithium heparin gel+, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 4,9 ml
Diameter 13 mm
Length including cap 106 mm
Length excluding cap 90 mm
Length including cap and plunger 127 mm

Material & colours

Colour of product transparent
Cap material High Density Polyethylene (HD-PE)
Colour of cap orange
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of inner packaging case
Piece(s) / inner box 50
Piece(s) / case 500
Piece(s) / pallet 36000
Depth of case 388 mm
Width of case 358 mm
Height of case 150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Lithium heparin gel+, 4.9 ml, cap orange, (LxØ): 90 x 13 mm, with paper label

Page 3

Case volume	0,0208 cbm
Weight of product	0,0057 kg
Weight of case	3,63 kg
EAN of inner box	4038917523158
EAN case	4038917523141

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT