



**Product Specification**  
**S-Monovette® Citrate 3.2%, 1.8 ml, cap green,  
(LxØ): 75 x 13 mm, with plastic label**

**Page 1**



### Product description

Order number	04.1955
Product description	S-Monovette® Citrate 3.2%, preparation: Trisodium citrate solution, 1.8 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label / print: transparent/green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Trisodium citrate solution
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/green
Graduation	fill mark
Cap	membrane screw cap
Type of cap	screw cap
Recommended centrifugation	3000 x g
Centrifugation temperature	20 °C
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-01-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**





**Product Specification**  
**S-Monovette® Citrate 3.2%, 1.8 ml, cap green,**  
**(LxØ): 75 x 13 mm, with plastic label**

**Page 3**

Case volume	0,0188 cbm
Weight of product	0,0048 kg
Weight of case	3,16 kg
EAN of inner box	4038917415880
EAN case	4038917415873

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-01-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**