



Product Specification

S-Monovette® neutral, 2.7 ml, cap white, (LxØ): 66 x 11 mm, with plastic label

Page 1

Product description

Order number	05.1729
Product description	S-Monovette® neutral, without preparation, 2.7 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 66 x 11 mm, with plastic label, label / print: black, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Type of cap	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20 °C
Centrifugation time	10 Min
Coding colour	white
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	2,7 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

Material & colours

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® neutral, 2.7 ml, cap white, (LxØ): 66
x 11 mm, with plastic label**

Page 2

Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 20000
Type of inner packaging case
Piece(s) / inner box 50
Piece(s) / case 500
Piece(s) / pallet 56000
Depth of case 333 mm
Width of case 303 mm
Height of case 125 mm
Case volume 0,0126 cbm
Weight of product 0,0038 kg
Weight of case 2,43 kg
EAN of inner box 4038917240895
EAN case 4038917240888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT