




Product Specification
S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with plastic label

Page 1



Product description

Order number	05.1729
Product description	S-Monovette® neutral Z, without preparation, 2.7 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 66 x 11 mm, with plastic label, label/print: white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Coding colour	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with plastic label**

Page 2

Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	2.7 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
66 x 11 mm, with plastic label

Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0038 kg
Weight of case	2.43 kg
EAN of inner box	4038917240895
EAN case	4038917240888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

