



Product description

Order number	51.5100.100
Product description	S-Monovette® POC Collect Kit, cap green, with paper label, label / print: green, set components: 1 S-Monovette® 1.2 ml Lithium Heparin (green), 1 POCT DISPENSER and 1 Multi-Adapter Luer Lock

Product characteristics

Label/ Print	with paper label
Colour of print/label	green

Material & colours

Colour of product	transparent
Colour of cap	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
Batched	yes

Packaging

Minimum order qty.	50
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	50
Piece(s) / pallet	4800
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0,0135 cbm
Weight of product	0,0091 kg
Weight of case	0,86 kg
EAN of inner box	4038917398763
EAN case	4038917398756

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

