



Product description

Order number	83.1850
Product description	Membrane pump, for pressure filtration, for Filtropur LS, 1 piece(s)/case

Product characteristics

Application	for pressure filtration
Suitable for	Filtropur LS

Size

Width of product	99 mm
Height of product	177 mm
Length of product	130 mm

Purity & certification

Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE supplier
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Membrane pump, for Filtropur LS

Page 2

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	50
Case volume	0.0157 cbm
Weight of product	3.159 kg
Weight of case	3.51 kg
EAN of inner box	4038917087940
EAN case	4038917006101

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT