



Product description

Order number	83.3952
Product description	Cell scraper, 2-position blade, size: L, handle length: 36 cm, blade length: 1.7 cm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Size

Length of product	36.9 cm
Length of blade	1.7 cm
Length of stem	36 cm
Width of blade	1.7 cm

Material & colours

Product material	Polystyrene (PS)
Colour of product	red/opaque
Stem material	Polystyrene (PS)
Blade material	Thermoplastic elastomers (TPE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	100
Piece(s) / pallet	7000
Depth of case	340 mm
Width of case	137 mm
Height of case	495 mm
Case volume	0.023 cbm
Weight of product	0.0164 kg
Weight of case	2.4 kg
EAN of inner box	4038917426695
EAN case	4038917426688

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

