



### **Product description**

Order number	83.3990.175
Product description	Screw cap, for tissue culture flask T-175, 2-position screw cap, without filter, red, TC Tested, 1 piece(s)/bag

### **Product characteristics**

Cap	2-position screw cap
Closure type	screw cap
Suitable for	tissue culture flask T-175
Filter	no

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	red
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, for tissue culture flask T-175, without filter**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	16800
Depth of case	288 mm
Width of case	218 mm
Height of case	150 mm
Case volume	0.0094 cbm
Weight of product	0.0054 kg
Weight of case	0.92 kg
EAN of inner box	4038917359986
EAN case	4038917359979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**