




Product Specification

**Serological pipette, with tip, plugged, 2 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 1



Product description

Order number	86.1252.001
Product description	Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Colour of print/label	green
Graduation	yes
Design	with tip, plugged
Volume of work	2 ml
Coding colour	green

Size

Diameter	6 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, with tip, plugged, 2 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Gamma irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging blister
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 54000
Depth of case 420 mm
Width of case 250 mm
Height of case 290 mm
Case volume 0.0323 cbm
Weight of product 0.0044 kg
Weight of case 5.03 kg
EAN of inner box 4038917127196
EAN case 4038917021081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

