




Product Specification
**Serological pipette, with tip, plugged, 25 ml,
sterile, 20 piece(s)/bag**

Page 1



Product description

Order number	86.1685.020
Product description	Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, 20 piece(s)/bag

Product characteristics

Colour of print/label	red
Graduation	yes
Design	with tip, plugged
Volume of work	25 ml
Coding colour	red

Size

Diameter	15 mm
Length of product	345 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, with tip, plugged, 25 ml,
sterile, 20 piece(s)/bag**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	20
Piece (s) / outer case	200
Piece(s) / pallet	10800
Depth of case	413 mm
Width of case	250 mm
Height of case	270 mm
Case volume	0.0279 cbm
Weight of product	0.0152 kg
Weight of case	3.61 kg
EAN of inner box	4038917133029
EAN case	4038917040785

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT