




Product Specification
quadriPERM® , square tissue culture dish, not TC-treated, 4 well, 12 piece(s)/bag

Page 1



Product description

Order number	94.6077.307
Product description	quadriPERM® , square tissue culture dish, not TC-treated, 4 well, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 12 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Number of wells	4
Growth area	24.9 cm ²
Volume of work	5 ml - 10 ml
Type of surface treatment	not TC-treated
Coding colour	green

Size

Width of product	89.5 mm
Height of product	17.7 mm
Length of product	130.7 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**quadriPERM® , square tissue culture dish, not
TC-treated, 4 well, 12 piece(s)/bag**

Page 2

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	48
Type of smallest subpackaging	bag
Piece(s) / inner box	48
Piece (s) / outer case	192
Piece(s) / pallet	3840
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0.0134 cbm
Weight of product	0.0366 kg
Weight of case	8.87 kg
EAN of inner box	4038917256131
EAN case	4038917256124

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-23

This document was prepared by EDP support and is valid without signature.



SARSTEDT