



### **Product description**

Order number	94.6077.434
Product description	flexiPERM® , con A, re-usable silicon insert, for dividing the growth surface, 1 well, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

### **Product characteristics**

Application	for dividing the growth surface
Number of wells	1
Growth area	1.1 cm <sup>2</sup>
Volume of work	1 ml - 1,5 ml

### **Size**

Diameter	12 mm
Height of product	8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	silicone
Colour of product	transparent

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	non-pyrogenic/endotoxin-free, non-cytotoxic
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	7200
Depth of box	185 mm
Width of box	100 mm
Height of box	30 mm
Depth of case	388 mm
Width of case	253 mm
Height of case	225 mm
Case volume	0.0221 cbm
Weight of product	0.0037 kg
Weight of case	1.42 kg
EAN of inner box	4038917255264
EAN case	4038917255257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

