



Product description

Order number	83.3922.400
Product description	Cell culture plate, 24 well, material: PS, surface: BIOFLOAT?, for spheroid culture, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Cell type	spheroid culture
Base shape	flat base
Number of wells	24
Growth area	1.82 cm ²
Volume of work	1 ml
Vessel type	Plate

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Max. volume	3.39 ml
Width of product	85.2 mm
Height of product	22.45 mm
Length of product	127.8 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Polystyrene (PS)

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	3
Type of smallest subpackaging	bag
Piece(s) / inner box	3
Piece (s) / outer case	18
Piece(s) / pallet	1080
Depth of box	205 mm
Width of box	140 mm
Height of box	80 mm
Depth of case	322 mm
Width of case	257 mm
Height of case	265 mm
Case volume	0.0219 cbm
Weight of product	0.0586 kg
Weight of case	2.16 kg
EAN of inner box	4066916042653
EAN case	4066916042660

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

