



**Product description**

Order number	65.386.007
Product description	Colour-coded inserts, for CryoPure tubes, violet, 100 piece(s)/bag

**Product characteristics**

Graduation	no
Suitable for	CryoPure tubes

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	violet

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	3000
Piece(s) / pallet	972000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0002 kg
Weight of case	0.68 kg
EAN of inner box	4038917106184
EAN case	4038917077729

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

