



Product Specification
CryoPure tubes, 1.2 ml, QuickSeal screw cap, white

Page 1



Product description

Order number	72.377
Product description	CryoPure tubes, 1.2 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, white, external thread, nominal volume: 1 ml, Cryo Performance Tested, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Closure type	screw cap
Internally or externally threaded cap	external thread
Volume of work	1 ml
Vessel type	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
CryoPure tubes, 1.2 ml, QuickSeal screw cap, white

Page 2

Size

Sample volume	1 ml
Max. volume	1.2 ml
Diameter	12.65 mm
Length of product	36.7 mm
Length excluding cap	29.8 mm

Material & colours

Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	104000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.0015 kg
Weight of case	4.24 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**CryoPure tubes, 1.2 ml, QuickSeal screw cap,
white**

Page 3

EAN of inner box	4038917081382
EAN case	4038917080507

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT