



Product description

| | |
|---------------------|--|
| Order number | 72.379.005 |
| Product description | CryoPure tubes, 2 ml, tube: PP, Quickseal screw cap, HD-PE, green, external thread, nominal volume: 1.8 ml, Cryo Performance Tested, 50 piece(s)/bag |

Product characteristics

| | |
|---------------------------------------|---------------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | yes |
| Cap | Quickseal screw cap |
| Type of cap | screw cap |
| Internally or externally threaded cap | external thread |
| Volume of work | 1,8 ml |
| Tube | Tube |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

| | |
|----------------------|----------|
| Sample volume | 1,8 ml |
| Max. volume | 2 ml |
| Diameter | 12,65 mm |
| Length of product | 46,7 mm |
| Length excluding cap | 39,8 mm |

Material & colours

| | |
|---------------|-----------------------------------|
| Cap material | High Density Polyethylene (HD-PE) |
| Colour of cap | green |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | IVD product |
| CE certified | CE |
| Purity standard | Cryo Performance Tested |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------|---------------|
| Minimum order qty. | 500 |
| Type of inner packaging | bag |
| Piece(s) / inner box | 500 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 104000 |
| Depth of box | 293 mm |
| Width of box | 183 mm |
| Height of box | 118 mm |
| Depth of case | 598 mm |
| Width of case | 378 mm |
| Height of case | 135 mm |
| Case volume | 0,0305 cbm |
| Weight of product | 0,002 kg |
| Weight of case | 5,17 kg |
| EAN of inner box | 4038917106054 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN case

4038917080583

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT