



### Product description

Order number	72.380.007
Product description	CryoPure tubes, 2 ml, tube: PP, Quickseal screw cap, HD-PE, violet, Internal thread, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	Quickseal screw cap
Type of cap	screw cap
Internally or externally threaded cap	Internal thread
Volume of work	1,6 ml
Tube	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Sample volume	1,6 ml
Max. volume	2 ml
Diameter	12,65 mm
Length of product	47,7 mm
Length excluding cap	37,8 mm

### Material & colours

Cap material	High Density Polyethylene (HD-PE)
Colour of cap	violet
O-ring material	silicon
Tube material	Polypropylene (PP)

### Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0,0305 cbm
Weight of product	0,0019 kg
Weight of case	4,99 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN of inner box                      4038917106108  
EAN case                                    4038917080644

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

