



### Product description

Order number	72.380.992
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, HD-PE, colour mix, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Closure type	screw cap
Internally or externally threaded cap	Internal thread
Design	with silicon O-ring
Volume of work	1,6 ml
Vessel type	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Sample volume	1.6 ml
Max. volume	2 ml
Diameter	12.65 mm
Length of product	47.7 mm
Length excluding cap	37.8 mm

### **Material & colours**

Closure material	High Density Polyethylene (HD-PE)
Colour of cap	colour mix
O-ring material	silicone
Tube material	Polypropylene (PP)

### **Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.0019 kg
Weight of case	4.96 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**CryoPure tubes, 2 ml, QuickSeal screw cap,  
colour mix**

EAN of inner box                      4038917228480  
EAN case                                 4038917228473

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

