



Product description

Order number	72.383.005
Product description	CryoPure tubes, 5 ml, tube: PP, QuickSeal screw cap, HD-PE, green, external thread, nominal volume: 4.5 ml, Cryo Performance Tested, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	QuickSeal screw cap
Closure type	screw cap
Internally or externally threaded cap	external thread
Volume of work	4,5 ml
Vessel type	Tube

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	4.5 ml
Max. volume	5 ml
Diameter	12.65 mm
Length of product	90.7 mm
Length excluding cap	83.8 mm

Material & colours

Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	250
Type of inner packaging	bag
Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	52000
Depth of box	293 mm
Width of box	183 mm
Height of box	118 mm
Depth of case	598 mm
Width of case	378 mm
Height of case	135 mm
Case volume	0.0305 cbm
Weight of product	0.0038 kg
Weight of case	5.05 kg
EAN of inner box	4038917106139

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN case

4038917080682

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT