



## Product description

Order number	93.873.481
Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, yellow, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 9 x 9, for 81 collection tubes, for CryoPure tubes 1.2?2.0 ml internal and external thread, 5 piece(s)/bag

## Product characteristics

Graduation	no
Cap	slip-on lid with ventilation function
Closure type	slip-on lid
Design	with numerical coding at each aperture
Suitable for	CryoPure tubes 1.2?2.0 ml internal and external thread
Storage capacity	81
Format	9 x 9

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Width of product	132 mm
Height of product	53 mm
Length of product	132 mm

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	yellow
Colour of cap	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	720
Depth of case	135 mm
Width of case	270 mm
Height of case	540 mm
Case volume	0.0197 cbm
Weight of product	0.1694 kg
Weight of case	3.91 kg
EAN of inner box	4038917122863
EAN case	4038917415095

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

