



## Product description

Order number	93.874.210
Product description	Cryobox, with numerical coding at each aperture, for low-temperature storage, material: PC, red, slip-on lid with ventilation function, cap: transparent, (LxWxH): 132 x 132 x 53 mm, format: 10 x 10, for 100 collection tubes, for CryoPure tubes 1.2 - 2.0 ml internal thread, 5 piece(s)/bag

## Product characteristics

Graduation	no
Cap	slip-on lid with ventilation function
Closure type	slip-on lid
Design	with numerical coding at each aperture
Suitable for	CryoPure tubes 1.2 - 2.0 ml internal thread
Storage capacity	100
Format	10 x 10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Size**

Width of product	132 mm
Height of product	53 mm
Length of product	132 mm

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	red
Colour of cap	transparent

### **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Batched	yes

### **Packaging**

Minimum order qty.	5
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	720
Depth of box	267 mm
Width of box	130 mm
Height of box	130 mm
Depth of case	550 mm
Width of case	280 mm
Height of case	145 mm
Case volume	0.0223 cbm
Weight of product	0.1593 kg
Weight of case	3.71 kg
EAN of inner box	4038917122894
EAN case	4038917415132

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

