




Product Specification
Tissue culture dish, (ØxH): 35 x 10 mm, surface:
Standard

Page 1



Product description

| | |
|---------------------|---|
| Order number | 83.3900 |
| Product description | Tissue culture dish, (ØxH): 35 x 10 mm, material: PS, surface: Standard, for Adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag |

Product characteristics

| | |
|----------------|---------------------|
| Closure type | slip-on lid |
| Cell type | Adherent cells |
| Base shape | SUREGrip Design |
| Growth area | 9.4 cm ² |
| Volume of work | 3 ml |
| Vessel type | Dish |
| Coding colour | red |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tissue culture dish, (ØxH): 35 x 10 mm, surface:
Standard

Page 2

Size

| | |
|-------------------|-------|
| Diameter | 35 mm |
| Height of product | 10 mm |

Material & colours

| | |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Closure material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | TC Tested |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 10 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 30000 |
| Depth of case | 322 mm |
| Width of case | 257 mm |
| Height of case | 245 mm |
| Case volume | 0.0203 cbm |
| Weight of product | 0.0038 kg |
| Weight of case | 2.59 kg |
| EAN of inner box | 4038917337953 |
| EAN case | 4038917365789 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

