



Product description

| | |
|---------------------|--|
| Order number | 83.3901.002 |
| Product description | Tissue culture dish, (ØxH): 60 x 15 mm, with grid, material: PS, surface: Standard, for Adherent cells, coding colour: red, TC Tested, 10 piece(s)/bag |

Product characteristics

| | |
|----------------|----------------------|
| Closure type | slip-on lid |
| Cell type | Adherent cells |
| Base shape | SUREGrip Design |
| Growth area | 22.1 cm ² |
| Volume of work | 5 ml |
| Vessel type | Dish |
| Coding colour | red |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|-------|
| Diameter | 60 mm |
| Height of product | 15 mm |

Material & colours

| | |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Closure material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | TC Tested |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 10 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 12000 |
| Depth of case | 332 mm |
| Width of case | 332 mm |
| Height of case | 380 mm |
| Case volume | 0.0419 cbm |
| Weight of product | 0.0059 kg |
| Weight of case | 4.02 kg |
| EAN of inner box | 4038917338103 |
| EAN case | 4038917338097 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

