



### Product description

Order number	83.3903.300
Product description	Tissue culture dish, (ØxH): 150 x 20 mm, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, TC Tested, 5 piece(s)/bag

### Product characteristics

Closure type	slip-on lid
Cell type	challenging adherent cells
Base shape	SUREGrip Design
Growth area	151.4 cm <sup>2</sup>
Volume of work	36 ml
Vessel type	Dish
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Diameter	150 mm
Height of product	20 mm

### **Material & colours**

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	1800
Depth of case	438 mm
Width of case	324 mm
Height of case	540 mm
Case volume	0.0766 cbm
Weight of product	0.0535 kg
Weight of case	6.92 kg
EAN of inner box	4038917338318
EAN case	4038917338301

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

