



Product Specification
**Cell culture flask, T-25, surface: Standard,
2-position screw cap**

Page 1



Product description

Order number	83.3910
Product description	Cell culture flask, T-25, material: PS, surface: Standard, for Adherent cells, coding colour: red, 2-position screw cap, TC Tested, 10 piece(s)/bag

Product characteristics

Cap	2-position screw cap
Cell type	Adherent cells
Growth area	24.7 cm ²
Volume of work	7 ml - 12,5 ml
Vessel type	flask
Coding colour	red

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Cell culture flask, T-25, surface: Standard,
2-position screw cap**

Page 2

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	300
Piece(s) / pallet	6000
Depth of case	595 mm
Width of case	408 mm
Height of case	350 mm
Case volume	0.085 cbm
Weight of product	0.0142 kg
Weight of case	6.61 kg
EAN of inner box	4038917337595
EAN case	4038917337588

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT