



Product description

Order number	83.3910.502
Product description	Cell culture flask, T-25, material: PS, surface: Suspension, for suspension cells, coding colour: green, Filter cap, TC Tested, 10 piece(s)/bag

Product characteristics

Cap	Filter cap
Cell type	suspension cells
Growth area	25 cm ²
Volume of work	7 ml - 12,5 ml
Vessel type	flask
Coding colour	green

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	green

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Cell culture flask, T-25, surface: Suspension,
Filter cap

Page 2

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	10
Piece (s) / outer case	300
Piece(s) / pallet	6000
Depth of case	602 mm
Width of case	408 mm
Height of case	340 mm
Case volume	0.0835 cbm
Weight of product	0.0172 kg
Weight of case	7.51 kg
EAN of inner box	4038917337809
EAN case	4038917337793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT