



**Product description**

Order number	83.3911.002
Product description	Cell culture flask, T-75, material: PS, surface: Standard, for Adherent cells, coding colour: red, Filter cap, TC Tested, 5 piece(s)/bag

**Product characteristics**

Cap	Filter cap
Cell type	Adherent cells
Growth area	75 cm²
Volume of work	21 ml - 55 ml
Vessel type	flask
Coding colour	red

**Material & colours**

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Colour of cap	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Purity & certification**

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	100
Piece(s) / pallet	1600
Depth of case	445 mm
Width of case	402 mm
Height of case	435 mm
Case volume	0.0778 cbm
Weight of product	0.05 kg
Weight of case	6.78 kg
EAN of inner box	4038917337687
EAN case	4038917366052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

